EN

FR

DE

ES

IT

NL

NO

JA

KO

Little Baby QCPR

Important Product Information

Laerdal

helping save lives

English

Little Baby Q CPR

This Important Product Information covers Little Baby Q CPR. Unless otherwise specified the information applies to all product configurations.

Warnings and Cautions

A Warning states a condition, hazard, or unsafe practice that can result in serious personal injury or death.

A Caution states a condition, hazard, or unsafe practice that can result in minor personal injury or damage to the product.

Notes

A note states important information about the product or its operation.

Care, Maintenance and Cleaning

Keep the manikin clean and in a hygienic condition. Replace the airway assembly after every class, especially if mouth to mouth resuscitation is practiced.

We recommend that you use a separate Laerdal Manikin Face for every student. If several students use one manikin face, thoroughly sanitize the manikin face after every use. Clean all skin parts regularly. Use warm soapy water or Laerdal manikin wipes. Use of Manikin Face Shields provides a clean barrier between your lips and the manikin face.

These can discolor the manikin:

- Pigments from lipstick and pens
- Latex gloves
- Using clothes other than what comes with the product originally.

Cautions

- *Do not use liquids to clean inside the manikin.*
- *Do not submerge any parts.*

Battery

Only use consumer grade, non-rechargeable, carbon zinc or alkaline batteries, LR6 1.5V AA. Batteries must be replaced with correct type. Remove batteries if product is not to be used for an extended period.

Symbol Glossary

| | |
|---|--|
|  | CE mark |
|  | Australian Radiocommunications and EMC Compliance Mark |
|  | KC (Korea Certification) certification mark |
|  | Japanese MIC Certification |
|  | WEEE mark |
|  | Manufacturer |
|  | Date of Manufacture |
|  | Reference Order Number |
|  | LOT Number |

Specifications

| | |
|-----------------------|---|
| Operating temperature | 5 °C to 35 °C (41 °F to 95 °F) Degradation in product performance may be experienced when used at low temperatures (< 15 °C / 59 °F) |
| Storage temperature | -10 °C to 50 °C (14 °F to 122 °F) |
| Humidity | 5 % - 95 % relative humidity |
| Weight | 1.8 kg (± 0.2 kg) |
| Dimensions | 65 x 26 x 14 cm (25 x 10 x 5 inches) |

Warranty

Refer to the Laerdal Global Warranty for terms and conditions. For more information visit www.laerdal.com

Federal Communications Commission Statement and Industry Canada Statements

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference
2. This device must accept any interference received, including interference that may cause undesired operation.

Caution: Changes or modifications not expressly approved by Laerdal Medical could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.
Increase the separation between the equipment and receiver:

Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

Consult the dealer or an experienced radio/TV technician for help.

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

1. This device may not cause interference.
2. This device must accept any interference, including interference that may cause undesired operation of the device.

Little Baby QCPR

Contains FCC ID: QHQ-QCPR-LB
IC: 20263-QCPRLB

Korea

R-C-Lm1-QCPR-LB

Japan

MIC certification 012-190010

EU

This product is in compliance with the essential requirements of Council Directive 2014/53/EU on Radio Equipment (RED), and Council Directive 2011/65/EU on Restriction on the use of certain hazardous substance (RoHS).

Waste Handling

The QCPR Sensor contains electronic components. Dispose of it at the applicable recycling facility in accordance with local regulations.

European directive 2012/19/EU (WEEE)
WEEE: this appliance is marked according to the European directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE). By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

The symbol on the product, or on the documents accompanying the product, indicates that this appliance may not be treated as household waste. Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment. Disposal must be carried out in accordance with local environmental regulations for waste disposal. For more detailed information about treatment, recovery and recycling of this product, please contact your local city office, your household waste disposal service or the Laerdal representative where you purchased the product.

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Reorient or relocate the receiving antenna.
Increase the separation between the equipment and receiver.
Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
Consult the dealer or an experienced radio/TV technician for help.

Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

1. This device may not cause interference.

2. This device must accept any interference, including interference that may cause undesired operation of the device.

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© 2019 Laerdal Medical AS. All rights reserved.

Manufactured in China for:
Laerdal Medical AS
P.O. Box 377
Tanke Svilandsgate 30, 4002 Stavanger, Norway
T: (+47) 51 51 17 00

Printed in China

20-15510 Rev B

www.laerdal.com

Laerdal
helping save lives