

The DefibSafe2 Lamppost Bracket



KEY FEATURES:

Greater Freedom of Placement:

With our bracket, our DS2 defibrillator cabinets can be conveniently placed on lampposts, expanding their accessibility and coverage to previously inaccessible areas.

Enhanced Visibility and Accessibility:

By mounting the DS2 Defibrillator cabinet on lampposts, they become highly visible, even during nighttime or adverse weather.

Optimal Space Utilisation:

Our bracket is designed to integrate with existing lamppost infrastructure making efficient use of urban space with additional clutter.

Easy Installation:

Our lamppost bracket is engineered for quick and hassle-free installation and is tested and approved for UMSUG (unmetered charging as used with streetlights).

Cost Effective:

Our bracket offers a cost-effective alternative to placing a cabinet on existing buildings and avoids the cost of maintenance work.

Community Engagement and Awareness:

Installing defibrillators on lampposts fosters a sense of community and well-being.

Public-Private Partnerships:

Local businesses and organizations can sponsor the installation of defibrillators on lampposts promoting their commitment to community welfare while improving public safety.

Compliance and safety:

Our bracket is manufactured to the relevant safety standards ensuring that it is a reliable solution for both rural and urban settings.

Empowerment and Peace of Mind:

By investing in our lamppost defibrillator bracket, communities can equip themselves with the tools to save lives.

ALSO AVAILABLE

iPAD SP1 Defibrillator
iPAD SPR Defibrillator
NFK200 Defibrillator



WEL
Medical
Life Matters

Contact us now
if you need more information or
if you would like to place an order:

T 01252 344007

orders@welmedical.com | www.welmedical.com

The DefibSafe2 Cabinet

Easy
to see

Keep your iPad SP1
secure until needed.

The DefibSafe2 ensures your defibrillator is easily located but remains safe.

Key features:

- LED lighting for high visibility at night
- Tough, weather proof casing (IP66 rated)
- Thermostatically controlled heater
- Unique location code for emergency services
- 10 year warranty (restrictions apply)
- Stainless steel lock (unlockable version available)
- Designed & Manufactured in the UK
- Tough Polyethylene body
- Polycarbonate Door
- Metal hinges with security fixings
- Double Skin Construction
- Low Energy – Low Voltage
- 3 Bolt fixing to wall
- Protection against the elements

Option for bespoke company logo on front of cabinet at an additional charge.

Lockable and non-lockable versions available.



DefibSafe 2 Cabinet

